



RNA – Scales & Yields

RNA Oligos			
	Scale(μmol)	OD-Yield (15-30mer)	OD-Yield (30-60mer)
RNA	0,04	> 2	> 1
	0,2	> 5	> 3
	1	> 10	> 8
2'-O-Methyl RNA	0,04	> 2	> 1
	0,2	> 5	> 3
	1	> 10	> 8
Duplex RNA	0,04	> 2	> 1
	0,2	> 5	> 3
	1	> 10	> 8

- for longer RNA please inquire
- Standard Purification Method: RNA'se free HPLC
- PAGE Purification on request
- All modifications are available – please see our overview under (Standards/Modifications)